



LONG LOOP SCISSORS

CODE

Right Hand with short round ended blade - LPR-1
Left Hand with short round ended blade - LPL-1

DESCRIPTION



Designed with an extended loop to give increased strength and control over the cutting action.

Assists children with poor motor control, particularly those with Dyspraxia.

HOW TO USE

Place the index finger outside the loop to provide guidance and control. The middle, ring and little fingers are inserted through the loop to provide the strength for the cutting action.

KEY FEATURES

- Stainless steel blades with guaranteed secure fixing to handles
- Extended loop allows strength and control of all fingers to be used
- Comfortable grip given by wide finger contact area
- Colour coded handles easily identify whether for right hand (blue), or for left hand (green)

RECOMMENDED FOR

Weak Grip	Poor Hand Control	Tremor



RELATED PRODUCTS

Also available as adult versions, please ask for details of LPR-5 and LPL-5
For self-opening versions ask for details of LPR-1/SO and LPL-1/SO

Peta (UK) Ltd

Mark's Hall, Mark's Hall Lane, Margaret Roding, Dunmow, CM6 1QT, England

Tel : +44 (0) 1245 231118 Fax : +44 (0) 1245 231811

E-mail : sales@peta-uk.com Internet : peta-uk.com